

For further assistance, please
visit [Tenovi.com/support](https://tenovi.com/support).

tenovi

© 2025 Tenovi Co.

All rights reserved. Republication, duplication or modification
in whole or in part is prohibited without prior written permission
from Tenovi. All trademarks are owned by or licensed to
Tenovi, Co. unless otherwise noted.

Revision 05.01.2024

1

Power Your Gateway

Plug your Gateway into a power outlet with the provided power cord and wait until the LED on the Gateway turns **RED**.



2


Pair Watch To Gateway

When your watch first turns on, select language. A screen will appear asking to **"Scan to Bind,"** just wait for the next screen to show up asking **"Agree to the phone binding request"** and **select the blue check mark.**



3 Charge your Tenovi Watch

Before first use, **plug the charging cable into a power supply with a USB-A port** and the other end to the watch which will snap on magnetically. **Recharge** every 5-7 days.




5 Wear The Tenovi Watch

Wear your watch every day unless instructed otherwise by your healthcare provider.


Do not pair your watch with any phone apps.

Do not submerge watch in water or wear in the shower. You may leave watch on while washing your hands.



4 Put On The Tenovi Watch

Place the band a finger width away from your wrist. **Adjust the strap tightness** so that it is not loose but still comfortable.



6 Transmit Data to your Provider

Every 12 hours, when you are within **25 feet** of the Gateway, the Tenovi Watch will **automatically connect** to the Tenovi Gateway and **transmit data** to your provider. The LED will turn **red** every morning and switch to **green** after the first connection in the day.



Measurement Needed



Measurement Successful