# Tenovi Blister Pack

## User Manual

tenovi

Model# : D-TE-BPCG-A

## Important Safety and Product Information:

### Warnings:

Failure to heed the warning below could result in an accident or medical event resulting in death or serious injury

This device contains electrical equipment that could cause injury if not handled properly.

Always consult your physician before beginning or modifying any healthcare program.

Do not allow children or pets to play with the device. The product contains small components that can be a choking hazard.

Do not attempt to replace any parts or open the enclosure. Doing so will void the warranty.

Do not use abrasive cleaners to clean your device.

### FCC Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference,

and (2) This device must accept any interference received, including interference that may cause undesired operation.

**Warning:** Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by removing the batteries from the equipment.

- Reorient or relocate the receiving antenna. Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## FCC Radiation Exposure Statement

This device has been evaluated to meet general RF exposure requirements. The device can be used in portable exposure conditions without restriction.

**FCC ID:** 2AZT4-TENPC001

## FCC Radiation Exposure Statement

This device complies with Industry Canada licence-exempt RSS standard(s).

Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device."

Le présent appareil est conforme aux CNR d'Indus trie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes :

(1) l'appareil nedoit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouil lage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement."

# Getting Started

## Step 1:  Plug in Gateway

A **Connect** the Gateway power cord to an electrical outlet.

B **Plug** the Gateway power cord into the back of the Gateway.

C The Gateway LED will first flash purple for up to ten minutes. It will then turn red indicating that it is successfully connected to a cellular network.

## Step 2: Take Your Medication

To **take your daily medication** as directed by your doctor.

To take your medications, **open the pill well** corresponding to the current date or time. Remove the medication and take it as instructed by your doctor.

## Step 3: Use Your Blister Pack Daily

The Blister Pack will automatically detect when you open a compartment. This information will be saved by your Blister Pack. Whenever you are within 30 feet of the Tenovi Gateway, the saved information from your Blister Pack will automatically be sent to your care provider. At this time, the cellular Gateway will beep and turn green indicating that data has been sent properly.

## Step 4: Transmit Open Events

The Tenovi Blister Pack will save all open events, even if you are away from the Tenovi Gateway. Whenever you open a pill compartment within 30 feet of the Tenovi Gateway, all stored open events will be sent to your care provider through the Tenovi Gateway. Your Gateway will turn from red to green and will beep if a medication event was sent that day. Your Gateway LED will turn red at midnight and then turn green and will beep once your first open event for the day has been successfully transmitted. If you open more pill wells in the same day they will be transmitted to the Gateway, but the color will not change.



**Warnings:**

Ensure that the space between the Blister Pack and the Gateway is free from obstructions and that they are in the same room when in use.

Tenovi Co.
1 Cate Street, Ste 100
Portsmouth, NH 03801
Email Contact: support@tenovi.com

tenovi